






SurgiTAIX AG

EndoTAIX A1

for Minimal Invasive Surgery


Automatic Endo-Camera Holder System...




...your precise and persistent assistance for MIS

... a research & development product...

... based on long lasting expertise in engineering ...





EndoTAIX A1

anthropomorphic modular design

20 N holding force (optional > 50 N)

working radius about 500 mm

fail-safe locking <20 msec

one-hand operation by using the button on the endo-camera or on the handle

sterile interface for the endoscope

fixation on the OR-table

draping with a conventional endo-camera drape

components optional autoclavable

VDE-approved safety & system concept

We would like to thank all clinical partners and the Helmholtz-Institute for BME of the Aachen University for their continuing interest and support in research and development

For further information please contact:



Pauwelsstraße 19
D-52074 Aachen, Germany
Tel: +49 (0)241 963-2145
E-mail: info@surgitaix.com
I-net: www.surgitaix.com

